



Laboratory Equipment Pty Ltd

1500 Degrees HIGH TEMPERATURE TUBE FURNACES



Features

- Silicon carbide elements
- Work tunes from 25mm to 90mm
- Double skinned stainless steel case and tubular frame
- Light weight ceramic fibre insulation
- Microprocessor based PID Controller
- Can be manufactured in a range of configurations
- Made to customer requirements

For Further Information contact Laboratory Equipment on

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au