



Laboratory Equipment Pty Ltd

1700 & 1800 Degrees HIGH TEMPERATURE HORIZONTAL TUBE FURNACES



Features

- **Molybdenum disilicide elements Alumina work tubes**
- **Double skinned stainless steel case and tubular steel frame**
- **Light weight ceramic fibre insulation**
- **Microprocessor based PID controller with over temperature alarm**
- **Can be manufactured in a range of configurations**
- **Made to customer requirements**

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au