



Laboratory Equipment Pty Ltd

1700 Degrees HIGH TEMPERATURE MUFFLE FURNACES



Features:

- Molybdenum disilicide elements
- Stainless steel double skinned construction
- Door face opens away from the operator
- Microprocessor based programmable digital controller and overtemperature alarm
- Phase angle thyristor power control and power transformer integrally mounted

Optional

- Gas outlets

Specifications:

Model	HT 04/17	HT08/17	HT16/17	HT32/17	HT64/17	HT117/17
Int. Dim. (mm) HxWxD	150X150X150	150X150X300	200X260X300	200X260X600	400X400X400	400X400X800
Ext. Dim. (mm) HxWxD	700X470X450	700X470X550	750X580X550	750X580X750	1050X1750X700	1050X1750X1100
Operating Temperature	1700°C					
Weight (kg)	150	200	200	300	600	700
Power(kW)	3.6	6	8	12	14	20

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au