



**Laboratory Equipment Pty Ltd**

## 1800 Degrees HIGH TEMPERATURE MUFFLE FURNACES



### Features

- Molybdenum disilicide elements
- Stainless steel double skinned construction
- Door face opens away from the operator
- Microprocessor based programmable digital controller and over temperature alarm
- Phase angle thyristor power control and power transformer integrally mounted
- Optional gas outlets

### Optional

- Gas outlets

### Specifications:

Model	HT05/18	HT07/18	HT10/18	HT 16/18
Int. Dim. (mm) WxDxH	150x150x200	150x150x300	200x200x250	200x260x300
Ext. Dim. (cm) WxDxH	720x470x550	720x470x600	720x520x570	720x580x620
Operating Temperature	1800°C			
Power (kW)	5.0	7.5	8.0	12
Weight (kg)	130	150	180	230

For Further Information contact Laboratory Equipment on: Email: [sales@labec.com.au](mailto:sales@labec.com.au)

Phone: 02 9560 2811 Fax: 02 9560 6131 [www.labec.com.au](http://www.labec.com.au)