

CALIBRATION/BOTTLE TESTING BATH (+100°C)

LABEC Calibration/Bottle Testing

Bath features stainless steel interior and exterior. Digital temperature controller. Up to 5 filling valves with individual flow levels. Automatic filling from hot water mains and powerful element provides for fast recovery where constant draining occurs. Circulated contents provides stable temperature, whilst gabled lid prevents loss of temperature.



SPECIFICATIONS

Model	LHF-40	LHF-80
Int. Dim. (cm) H x L x W	45 x 62 x 37	56 x 70 x 45
Temp. Min	Ambient +5°C	
Temp. Max	100°C	
Storage Capacity (L)	40	80
Power (W)	3600	4000
Weight (kg)	45	52

Serving the Scientific, Medical and Research industries since 1945

LABEC



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd
email: sales@labec.com.au
Ph: 02 9560 2811 • Fax: 02 9560 6131
www.labec.com.au