

CHEMICAL STORAGE CABINET

LABEC Chemical Storage Cabinet

are purpose built for the storage of corrosive substances. The cabinets are manufactured to AS3780 and incorporate the following features:

- * Self closing door.
- * 150mm deep liquid tight sump.
- * Sump capable of containing 25% of storage capacity.

Please note:

These cabinets are designed for the safe storage of corrosives in sealed and closed containers.

Good housekeeping practices are essential to ensure a long service life.



Corrosive substances cabinets are manufactured to AS3780 and incorporate the FOLLOWING FEATURES:

- Self closing door.
- 150mm deep liquid tight sump.
- Sump capable of containing 25% of storage capacity.

SPECIFICATIONS

Model	P-HA3011	P-HA3016	P-HA3017	P-HA3012	P-HA3018
Int. Dim. (mm) H x W x D	520 x 835 x 530	520 x 420 x 375	785 x 420 x 375	1015 x 1015 x 375	1545 x 1015 x 420
Ext. Dim. (mm) H x W x D	800 x 920 x 615	800 x 505 x 460	1065 x 505 x 460	1295 x 1100 x 460	1825 x 1100 x 500
Storage Capacity	100 litres/kgs	30 litres/kgs	60 litres / kgs	160 litres/kgs	250 litres/kgs
Shelves	1	1	2	2	3

Please note:

These cabinets are designed for the safe storage of corrosives in sealed and closed containers. Good housekeeping practices are essential to ensure a long service life.

* AS3780 specifies a maximum of 50kg or L for packaging group 1 corrosive substances

Serving the Scientific, Medical and Research industries since 1945

LABEC



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au