

EGG INCUBATOR (+50°C)

LABEC'S Egg Incubator can control temperature, humidity, ventilation and egg cycle time.

With the capacity to store up to 88 eggs, the Labec Egg Incubator is fully adjustable to meet your requirements.



SPECIFICATIONS

Model	X-XYA-1
External Dimension (cm)	70 x 55 x 80
Capacity (L)	Up to 88 eggs
Power (W)	200W
Temp. Display Range	5-50°C
Measurement Accuracy	+/-0.1°C
Temp. Control	+/-0.1°C
Humidity Display Range	0-99% RH
Humidity Display Control	+/-5% RH
Over Egg Cycle	0.1-99.9hr (default settings 1.5hrs)
Eggs Over Time	1-255 seconds adjustable (default setting 180 seconds)
Ventilation Cycle	5-999 minutes (adjustable)
Ventilation Time	0-999 seconds (adjustable)
Relative Humidity	Less than 85%

Serving the Scientific, Medical and Research industries since 1945

LABEC



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au