

FLOTATION BATH (+65°C)

LABEC Flotation Bath is manufactured in Australia with high quality stainless steel for easy maintenance, excellent durability. With three models available there is a floatation bath for all tissue preparation. Available in either a round or square lift out black aluminum tray or the FB/L1 is an illuminated model available with glass tray.

The illuminated model has a high temperature lift out glass tray and a rear mounted 9 watt high efficiency lamp.



FEATURES:

- All stainless steel construction
- Digital temperature controllers
- Built-in side warming plate
- Round or Square anodised aluminium trays
- FB/L1 is an illuminated model available with glass tray (as shown above)

SPECIFICATIONS

Model	FB1	FB2	FB/L1
Int. Dim. (cm) H x W x D	25 x 20 x 7.5	20dia x 7.5	21 x 21 x 6
Ext. Dim. (cm) H x W x D	15 x 45 x 31	15 x 45 x 31	15 x 45 x 31
Temp. Min.	Ambient +5°C		
Temp. Max	65°C		
Capacity (L)	4	4	0.5
Power (W)	250		
Weight (kg)	10	10	10

Serving the Scientific, Medical and Research industries since 1945

LABEC



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au