

IMMERSION CIRCULATOR (+100°C)

LABEC Immersion Circulator

available with either digital set and display of actual operating temperatures or analogue settings. The PID controller provides very precise temperature control and is able to adapt to many different containers. Equipped with additional safety cut-off. Options Pump set for external temperature applications and S - Optional Safety Cut Off.



Model TEP-4 shown above

SPECIFICATIONS

Model	M-TEP-1	M-TEP-4
Ext. Dim. (cm) W x H x D	20 x 24 x 9	
Temp. Min	Ambient +5°C	
Temp. Max	100°C	
Pumping Capacity	12 litres per minute	
Power (W)	1000	
Immersion Depth Minimum	6cm	
Immersion Depth Maximum	15.5cm	
Controller	Analogue	Digital
Weight (kg)	2.3kgs	

Serving the Scientific, Medical and Research industries since 1945

LABEC



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd
email: sales@labec.com.au
Ph: 02 9560 2811 • Fax: 02 9560 6131
www.labec.com.au