

PHARMACY VACCINE REFRIGERATOR (+2°C to +8°C)

LABEC offers a range of glass door and solid door **Vaccine/Pharmacy refrigerators** which are available in smaller under bench models to larger units up to 440 litres capacity. The Labec Pharmacy Vaccine refrigerators are used for the purpose of storing vaccine products safely.

Approved by the Quality Care Pharmacy Program and The Pharmacy Guild of Australia the refrigerators have a temperature range of +2°C to +8°C with a temperature fluctuation of +/-2°C. The temperature can be precisely controlled using a PID Microprocessor controller which has Hi/Low audio and visual alarms and can be optioned for BMS monitoring.

FEATURES INCLUDE:

- Fully compliant with Quality Care Pharmacy Program (QCCP) Cold Chain Standards.
- External Digital temperature controller.
- CFC free refrigerant.
- Audible/Visual/BMS signals.
- Hi/Lo Temperature memory.
- Hi/Lo Temperature alarm.
- External alarm facility for connection to building maintenance system.
- Adjustable/removable shelves.
- Internal light.
- Keylock door.
- Tempered thermoglass door on glass door model.



APPROVED BY THE
Quality Care Pharmacy Program
and The Pharmacy Guild of Australia.



Serving the Scientific, Medical and Research industries since 1945

LABEC

LABEC

Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au

SPECIFICATIONS

PHARMACY VACCINE REFRIGERATOR (+2°C to +8°C)

Model	PVR-130	PVR-130GD	PVR-374	PVR-440*
Capacity	130	130	374	440
Int. Dim. (cm) H x W x D	70x44x44	66x46x46	150x43x50	170 x 53x 50
Ext. Dim. (cm) H x W x D	88x54x54	85x60x60	193x60x65	210x68x65
Temp. Controller	Digital Controller			
Temp. Min	+2°C			
Temp. Max	+8°C			
Temp. Setting Accuracy	+/- 0.1°C			
Temp. Fluctuation	+/- 2°C			
Refrigerant	CFC Free			
Interior Material	Eeasy Clean Plastic Formed Lining			
Exterior Material	Powder Coated Steel			
Door Material	Solid	Glass	Solid or Glass	Glass
Shelves (adjustable)	3	3	4	5
Options	Cable entry port, extra shelves, low temperature safety thermostat.			
Power (W)	200		220	300
Power Supply	240V / 50Hz			

* PVR440 pending approval with QCCP

Serving the Scientific, Medical and Research industries since 1945

LABEC



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au