


TISSUE FLOAT BATH (+60°C)

LABEC Tissue Float Bath is manufactured with high quality aluminum, used for applications such as stretching and drying cut tissues for histological, pathological, chemical, clinical and bacteriological laboratories.

A digital controller ensures precision temperature controller within the temperature range of Ambient +5°C to 60°C and the temperature is set by using soft touch key pads.

SPECIFICATIONS

| | |
|-------------------|--------------|
| Model | J-NBT |
| Size Dia x H (mm) | 220x60 |
| Temp. Min | Ambient +5°C |
| Temp. Max | 60°C |
| Heater(sheath) W | 300 |
| Weight (kg) | 2.75 |



Serving the Scientific, Medical and Research industries since 1945

LABEC

LABEC

Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd
email: sales@labec.com.au
Ph: 02 9560 2811 • Fax: 02 9560 6131
www.labec.com.au