



















Laboratory Equipment Pty Ltd

Ultra Low Temperature Freezer Accessories

Model:	DF8510	DF8514	DF8517	DF8520	DF8524
Amount of racks: (view from afront)	 12 racks	 16 racks	 16 racks	 20 racks	 20 racks
Type of racks for 2 inch boxes:	 3 deep, 5 high Part: UR0066	 3 deep, 5 high Part: UR0066	 4 deep, 5 high Part : UR0068	 4 deep, 5 high Part : UR0068	 5 deep, 5 high Part : UR0070
Type of racks for 3 inch boxes:	 3 deep, 3 high Part: UR0067	 3 deep, 3 high Part: UR0067	 4 deep, 3 high Part : UR0069	 4 deep, 3 high Part : UR0069	 5 deep, 3 high Part : UR0071
Max. of 2 inch:	180 boxes	240 boxes	320 boxes	400 boxes	500 boxes
Max. of 3 inch:	108 boxes	144 boxes	192 boxes	240 boxes	300 boxes



External Rack Dimensions:
(mm D x W x H)

425x145x283	425x145x283	558x145x283	558x145x283	692x145x283
-------------	-------------	-------------	-------------	-------------




Internal Rack Dimensions
(mm D x W x H)

For 2 inch boxes, each opening: 145 x 136 x 51

For 3 inch boxes, each opening: 145 x 136 x 88

6.1 Storage Racks for Upright Freezer

6.2 Storage Racks for Chest Freezer

Model:	DF9007	DF9010	DF9014	DF9017	DF9020
Amount of racks: (view from above)	 12 racks	 18 racks	 24 racks	 30 racks	 36 racks
Type of racks for 2 inch boxes:	 12 high Part: CR0062	 12 high Part: CR0062	 12 high Part: CR0062	 12 high Part: CR0062	 12 high Part: CR0062
Type of racks for 3 inch boxes:	 8 high Part: CR0063	 8 high Part: CR0063	 8 high Part: CR0063	 8 high Part: CR0063	 8 high Part: CR0063
Max. of 2 inch:	144 boxes	216 boxes	288 boxes	360 boxes	432 boxes
Max. of 3 inch:	96 boxes	144 boxes	192 boxes	240 boxes	288 boxes



External Rack Dimensions:
(mm D x W x H)

146x146x660	146x146x660	146x146x660	146x146x660	146x146x660
-------------	-------------	-------------	-------------	-------------

Internal Rack Dimensions
(mm D x W x H)

For 2 inch boxes, each opening: 145 x 138 x 51

For 3 inch boxes, each opening: 137 x 134 x 78

6.3 Fireboard Storage Boxes and Dividers



Thick Plasticized Card box Material

Easy to remove cover

External dimensions:

External dimensions of 2 inch boxes (model PP0060): 133x133x49 mm (DxWxH)

External dimensions of 3 inch boxes (model PP0061): 133x133x77 mm (DxWxH)

Internal Dimensions:

Internal dimensions of 2 inch boxes (model PP0060): 124x124x45 mm (DxWxH)

Internal dimensions of 3 inch boxes (model PP0060): 124x124x70 mm (DxWxH)

Possible dividers:

Dividers are also made of plasticized card box. Height of all dividers is 35 mm

Model: Description:

PP1049 Storage box divider 7 x 7

PP1064 Storage box divider 8 x 8

PP1081 Storage box divider 9 x 9

PP1100 Storage box divider 10 x 10

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131