

VORTEX TUBE MIXER

LABEC'S Vortex Tube Mixer is Compact, robust, & durable design for efficiency, with variable speed control up to 3000rpms, continuous operation or automatic press to start. A fast, thorough and quiet vortex mixing action, but vibration free.

Using high quality motor to achieve the performance the VM1000 has a single cup head suitable for mixing the content of one tube.



SPECIFICATIONS

Model	M-VM1000
External Dimension (mm) WxDxH	120x185x160
Speed Range	3,000 rpms (max)
Motor Hp	1/15HP
Material	Fibre glass reinforced plastic
Function	Speed Controller and timer
Capacity (L)	Rubber single cup head with Ø26 mm
Weight (kg)	3

Serving the Scientific, Medical and Research industries since 1945

LABEC



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au