

SPECIFICATIONS

AUTOCLAVE - BENCH MODEL

B & S CLASS - HOLLOW & SOLID INSTRUMENTS

Model	MES-8BS	MES-12BS	MES-17BS	MES-22BS
Int. Dim. (cm) Dia x Length	17 x 31	20 x 36	25 x 35.5	25 x 45
Ext. Dim. (cm) L x W x D	56 x 41 x 40.5	56 x 41x 40.5	59 x 46.5 x 45	59 x 46.5 x 45
Temp.range	105° to 134°C			
Capacity (L)	8	12	17	22
Power (W)	1800	1800	1800	2300
Power Source	10 AMP/ 240 Volt			
Controller	LCD Display and 7 step programmable controller			
Printer	Optional	Optional	Standard	Standard
Safety Features	Automatic door lock under pressure, safety valve, temperature control and pressure protection system together with a temperature protector mounted on the chamber wall.			
Weight (kg)	39	42	48	50

Serving the Scientific, Medical and Research industries since 1945

LABEC

LABEC

Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au