



Laboratory Equipment Pty Ltd

BENCH TOP SHAKERS



Used for culturing microbes and further culturing, mixing, extracting and manipulating microbiological and plant tissues. These benchtop shakers are designed for rotation type shaking, slide shaking, reciprocating and S movement.

Options

- Flask clamps 100ml, 250ml,500ml
- Flask clamps spring holders
- Spring rackslide rack
- spring rack
- universal rack

Specifications

Model	J-USR	J-NSS	J-TUS
Overall Dim. (mm) WxDxH)	300x350x120	300x350x120	450x380x150
Platform Size (mm) WxD	300x300	300x300	400x300
Action Type	rotation type	slide type	reciprocating,& rotation movement
Shaking Speed	35-250 rpm		
Capacity	16 x 100ml /9 x 250ml /6 x 500ml		15 x 250ml
Weight (kg)	23	21	20.3

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au