



Laboratory Equipment Pty Ltd

CALIBRATION/ BOTTLE TESTING BATH



Featuring stainless steel interior and exterior. Digital temperature controller. Up to 5 filling valves with individual flow levels. Automatic filling from hot water mains and powerful element provides for fast recovery where constant draining occurs. Circulated contents provide stable temperature, whilst gabled lid prevents loss of temperature.

Specifications

Model	LHF-40	LHF-80
Ext. Dim. (cm) HxLxW	45x62x37	56x70x45
Temp. Min	Ambient+5°C	
Temp. Max	100°C	
Storage Capacity (L)	40	80
Power (W)	3600	4000
Weight (kg)	45	52

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au