



**Laboratory Equipment Pty Ltd**

## Ceramic Laboratory Crucibles

### Standard Crucibles

(D=top diameter, L=length, d=bottom diameter, W=width, H=height, in mm)

Catalogue No.	Type	D (L)	d (W)	H	capacity (cm <sup>3</sup> )
---------------	------	-------	-------	---	-----------------------------

#### Combustion Boats:

CB-0001		(48)	(11)	7	1
CB-0005		(90)	(15)	10	5
CB-0025		(100)	(24)	15	25
CB-0200	Half cylinder	(146)	(70)	35	200

#### Conical Crucibles:

CC-0007	Conical	30	20	20	7
CC-0008	Conical	26	16	29	8
CC-0030	Conical	36	28	43	30
CL-0030	Low Form	48	20	30	30
CT-0040	Tall Form	48	25	46	40

**Cylindrical Crucibles:**

CN-0002		14		25	2
CN-0004-W		18		25	4
CN-0005-W		20		30	5
CN-0007		18		40	7
CN-0025		30		45	25
CN-0050		35		65	50
CN-0090		40		90	90
CN-0175		64		72	175
CN-0200		52		130	200
CN-0325		65		140	325
CN-0500		70		150	500
CN-0550		83		100	550
CN-0750		83		160	750
CN-1650		115		180	1650
CN-4000		170		210	4000

**Crucible lids**

CLS	flat disc	23		3	
CLM	inverted dish	45	40	4	
CLL	" "	91	85	9	

**Crystal Growing Crucibles (30 deg):**

CG-0020		24		100	20
CG-0090		36		150	90

**Circular Dishes:**

CD-0020		45		20	20
CD-0050		60		20	50
CD-0070		102		14	70
CD-0080		78		30	80
CD-0325		102		49	325
CD-1000		180		56	1000

**Hermetic Enclosures:**

HE-0250		85	62	89	250
---------	--	----	----	----	-----

**Rectangular Trays:**

RT-0050		(82)	(55)	17	50
RT-0200		(140)	(55)	44	200
RT-0700		(280)	190)	20	700

For Further Information contact Laboratory Equipment on: Email: [sales@labec.com.au](mailto:sales@labec.com.au) Phone: 02 9560 2811 Fax: 02 9560 6131 [www.labec.com.au](http://www.labec.com.au)