



Features

A constant temperature water bath fitted with circulator pump for more even temperature uniformity.

For Stirred applications, or where uniformity of temperature is imperative.

All stainless steel construction

Digital PID solid state temperature controllers

Standard models from 8 to 44 litres

- Custom sizes also available up to 170 Litres

Specifications

Model	NBCT2S	NBCT5S	NBCT7S	NBCT9S	NBCT11S
Capacity (litres)	8	12	16	25	44
Int. Dim. (cm) LxWxD	30X23X12	30X23X17	32X30X17	50X30X17	59X49X17
Ext. Dim. (cm) LxWxD	40X37X27	40X37X27	40X38X27	67X38X27	75X52X27
Temp. Controller	Digital PID Controller				
Temp. Min	Ambient +5°C				
Temp. Max	100°C				
Temp. Setting Accuracy	0.1°C				
Temp. Fluctuation	+/- 0.2°C				
Temp. Sensor	Pt 100ohms				
Heating System	Electric Element				
Interior Material	Stainless Steel				
Exterior Material	Stainless Steel				
Options	Flat or gabled lift of lids in either stainless steel or polycarbonate, manual reset over heat cut off, top mounted stirrer, large capacity models, refrigerated models.				
Safety Features	Over temperature alarm, temperature protection cut off, earth leakage breaker (ELB).				
Power (W)	750	750	1200	2000	2400
Weight (kg)	10	15	20	22	33