



Laboratory Equipment Pty Ltd

FLOTATION BATH



Features

- All stainless steel construction
- Digital temperature controllers
- Built-in side warming plate
- Round or Square anodised aluminium trays
- FB/L1 is an illuminated model available with glass tray (as shown above)

Specifications

Model	FB1	FB2	FB/L1
Int. Dim. (cm) HxWxD	25x20x7.5	25diax7.5	21x21x6
Ext. Dim. (cm) HxWxD	15x45x31	15x45x31	15x45x31
Temp. Min	Ambient+5°C		
Temp. Max	65°C		
Capacity (L)	4	4	0.5
Power (W)	250		
Weight (kg)	10	10	10

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au