

LABORATORY REFRIGERATOR

(+2°C to +8°C)

LABEC stocks 4 models of **Economy Laboratory Refrigerators** ranging from 70 litres to fit under the bench up to 370 litres. Labec laboratory refrigerators are designed for general purpose storage of non-critical products where alarms and monitoring are not generally required. These units are set up as manual defrost.

All use CFC Free urethane insulation bonded between painted steel exterior and white moulded plastic "easy clean" interior with curved corners. All models feature self closing solid doors, adjustable open wire white plastic coated full width shelves which are easily removable. A hermetically sealed manual defrost refrigeration system utilizing CFC-Free refrigerant provides. Power requirements for all models are 10A 240V AC 50Hz power supply.

If the bench height is not standard, the GLR-130 can be customised to a height range of 845mm-880mm at no extra cost.

* Manual defrost requires regular defrosting which must be carried out to maintain temperature uniformity. Long term storage is not possible with manual defrosting, defrosting may be required as often as monthly depending on environmental variables and sample loads. Auto defrosting available on all models upon request.



Serving the Scientific, Medical and Research industries since 1945

LABEC



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au

SPECIFICATIONS

LABORATORY REFRIGERATOR (2°C to 8°C)

Model	GLR-70	GLR-130	GLR-240	GLR-370
Int. Dim. (mm) H x WxD	500 x 380 x 350	700 x 440 x 400	1280 x 450 x 400	1765 x 545 x 600
Ext. Dim. (cm) H x W x D	800 x 440 x 470	880 x 540 x 540	1500 x 550 x 540	1825 x 595 x 640
Temp. Min	+2°C			
Temp. Max	+8°C			
Capacity (L)	70	130	240	370
Number of Shelves	3	3	4	5
Temperature Display	Digital			
Door	Solid	Solid	Solid	Solid or Glass
Power (W)	200	220	250	350
Weight (kg)	20	24	39	85

Serving the Scientific, Medical and Research industries since 1945



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au