

GENERAL PURPOSE OVEN

FAN FORCED (+300°C)

LABEC General Purpose Fan Forced Oven

is manufactured in Australia with high quality stainless steel for easy maintenance and excellent durability. Suitable for applications where general warming or heating is required.

Using fan forced convection heating which circulates the warm air inside the chamber thus providing improved uniformity and fast recovery after the door is opened.

A digital microprocessor based PID temperature controller controls the temperature within the range ambient +5°C to +300°C, temperature fluctuation at 105°C is +/-1°C, when optioned the controller can be connected via RS232/RS485 communication ports which can be connected directly to a data logger or computer to collect this information. Software is available for programming and monitoring.

To suit your functionality the addition of options such as a safety toughened glass viewing window, programmable controllers, timers and cable access ports are also available.



RIGHT: Model shown with glass viewing window option.

OTHER STANDARD FEATURES INCLUDE:

- RCD/ ELB (earth leakage breaker)
- Over Temperature Thermostat
- Stainless Steel exterior
- Type K Thermocouple
- Adjustable/removable shelving



SPECIFICATIONS

GENERAL PURPOSE OVEN

FAN FORCED (+300°C)

Model	ODWFH14	ODWFH18	ODWFH24	ODWFH36
Capacity (Litres)	32	66	160	300
Int. Dim. (cm) H x W x D	35 x 30 x 30	46 x 38 x 38	60 x 50 x 50	90 x 60 x 50
Ext. Dim. (cm) H x W x D	65 x 38 x 48	70 x 47 x 56	90 x 60 x 69	125 x 72 x 70
Temp. Min	Ambient+5°C			
Temp. Max	300°C			
Temp. Setting Accuracy	+/-0.1°C			
Temp. Fluctuation	at 105°C is +/- 1°C			
Temp. Sensor	Type K			
Heating System	Electric Element			
Interior Material	Stainless Steel			
Exterior Material	Stainless Steel			
Shelves	2			
Options	Glass or acrylic viewing window, castors, cable entry port, timer, programmable temperature profiler, additional shelves			
Safety Features	Over temperature alarm, temperature protection cut off, earth leakage breaker.			
Power Supply	240V/50Hz			
Power (W)	2400		3600	4800
Weight (kg)	28	38	48	94

Serving the Scientific, Medical and Research industries since 1945



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au