

# HEATING MANTLE (+450°C)

**LABEC** high quality **Heating Mantles** are available in 50ml, 100ml, 250ml, 500ml, 1 to 5 litres and 10 litres capacity.

The heating mantle has a built in simmerstat controller to regulate temperature to 450°C.

The flexible coiled heating element provides excellent heat transfer while absorbing shock to minimise flask damage.

The rigid housing supports the weight of the vessel and provides uniform heating of the vessel without the danger of an open frame.



## SPECIFICATIONS

Model - Manual Control	Capacity	Power
M-MN-50C	50ml	60W
M-MN-100C	100ml	80W
M-MN-250C	250ml	200W
M-MN-500C	500ml	300W
M-MN-1000C	1000ml	400W
M-MN-2000C	2000ml	500W
M-MN-3000C	3000ml	600W
M-MN-5000C	5000ml	700W
M-MN-10LC	10Litres	1300W

Serving the Scientific, Medical and Research industries since 1945

**LABEC**

**LABEC**

Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: [sales@labec.com.au](mailto:sales@labec.com.au)

Ph: 02 9560 2811 • Fax: 02 9560 6131

[www.labec.com.au](http://www.labec.com.au)