



Laboratory Equipment Pty Ltd

HEATING MANTLE



A range of heating mantles to suit flask type are available in 50ml, 100ml, 250ml, 500ml, 1 to 5 litres and 10 litres capacity.

Features

- Maximum operating temperature of 450°C
- Rigid aluminium casting
- Fabric mantles for effective heating and protection of glass vessels from thermal shock
- Manual voltage control from 0-100%

Specifications

Model- Manual Control	Capacity	Power
M-MN-50C	50ml	60W
M-MN-100C	100ml	80W
M-MN-250C	250ml	200W
M-MN-500C	500ml	300W
M-MN-1000C	1000ml	400W
M-MN-2000C	2000ml	500W
M-MN-3000C	3000ml	600W
M-MN-5000C	5000ml	700W
M-MN-10LC	10Litres	1300W

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au