



Features

Our high quality magnetic stirrer allows heating and stirring to be performed either separately or simultaneously.

The top plate is aluminium die casting coated with ceramic and is very resistant to corrosive chemicals.

Accurate temperature control with temperature thermostat.

The hot plate is separated from the main body by a heat resistant plate to protect the electronics and the controllers.

Specifications

Model	J-NMS
Stirring Speed	100-1600/min
Stirring Capacity	100 ml - 3 litres
Stirring Control	Electronic
Plate Area (mm)	180 x 180
Overall Dimension (cm) HxWxD	120x200x225
Hot Plate Temp. Max	up to 400°C
RPM (Maximum)	1600 per minute
Wattage (W)	630
Weight (kg)	3.65kgs