



**Laboratory Equipment Pty Ltd**

## HOT PLATE



**Our quality hot plates provide even stable heat using a ceramic coated top plate, coated with white for easy observation during use and thermostat for temperature control.**

### Specifications

Model	J-NHP
Plate Dimensions (mm) DiaxDPTH	180mm x 180mm
Overall Dimension (cm) HxWxD	120x200x225
Hot Plate Temp. Max	up to 400°C
Wattage(W)	600
Weight (kg)	2.5

For Further Information contact Laboratory Equipment on:

Email: [sales@labec.com.au](mailto:sales@labec.com.au)

Phone: 02 9560 2811

Fax: 02 9560 6131

[www.labec.com.au](http://www.labec.com.au)