



Laboratory Equipment Pty Ltd

ICE MACHINE



Flake type ice maker using a unique freezer-flake mechanism to provide clean, hard bits of ice. An automatic system from water supply to ice production. Ideal for medical and laboratory facilities.

Specifications

Model	IF300-150
Ext. Dim. (mm) WxDxH	700x700x975
Type	Flake
Ice Production	120kgs/day
Ice storage	90kgs
Refrigeration system	1/3Hp Compressor
Condenser	Air Cooled Condensing Unit
Power Requirements	240V / 10AMPS
Cabinet material	High Density polyethylene
Weight (kg)	114kgs

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au