



Laboratory Equipment Pty Ltd

LARGE CAPACITY FURNACES



Features:

- **Stainless steel construction**
- **Designed for to suit applications such as mechanical testing apparatus**
- **Up to 1800°C**
- **Rapid heat up time**
- **Light weight ceramic fibre insulation for rapid thermal cycling**
- **Digital microprocessor based PID temperature controller**
- **All units custom made to suit your requirements**
- **Optional multi-stage fully programmable controllers available with RS232/485**

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au