



Laboratory Equipment Pty Ltd

LARGE CAPACITY - DRUM HEATING OVEN



- Reinforced internal flooring and supports to suit 205 litre drum units.
- Suitable for applications where high levels of moisture are present.
- Stainless steel interior fully sealed to prevent leakage into the insulation layer
- Stainless steel exterior
- Full insulation on all surfaces to keep exterior cool.
- 7 day on/off timer
- Manual reset thermostat
- Top mounted 50mm extraction vents (for connection to flue system if required)
- Continuous positive ventilation through the top vents
- Alarm over temperature cut-off protection
- Digital microprocessor based PID temperature controller accuracy 0.1°C
- Temperature stability +/-1°C
- Options include higher operating temperature, chart recorder and programmable temperature controllers, stainless flue system to suit your building.

Specifications

Model	ODDF-2x205	ODDF-4x205
Int. Dim. (cm) HxWxD	100x135x70	100x135x135
Ext. Dim. (cm) HxWxD	130x155x85	130x170x170
Temp. Max	150°C*	150°C*
Capacity (Drums)	2x 205 litre drums	4x 205 litre drums
Capacity (L)	945	2400
Controller	PID Controller	PID Controller
Power (W)	3600	9000
Weight (kg)	600kgs	900kgs

* Optional temperature maximum 250°C

For Further Information contact Laboratory Equipment on: Email: sales@labec.com.au Phone: 02 9560 2811 Fax: 02 9560 6131 www.labec.com.au