



Laboratory Equipment Pty Ltd

MULTI PURPOSE WATER BATH



This water bath has versatile usability and applications. Purchased as a full set or as a basic water bath with optional attachments. By attaching or detaching the modular functional parts such as immersion circulating bath, Soxhlet water bath and shaking water bath to the basic unit, this can be used for multi purposes.

Features

- Shaking speed control system (digital setting/display with feedback system)
- Indirect heating system

Options

Programmable control system. This control system enables continuous operation and programmed operation where the temperature can be increased, maintained or decreased according to requirements.

Specifications

Model	J-MPB2
Bath Inner Size(mm) HxWxD	220x580x380
Overall Size (mm) HxWxD	440x820x540
Shaking Type	Reciprocating type
Temp. Min	Ambient+5°C
Temp. Max	95°C
Temp .Accuracy	+/- 0.1°C at 50°C
Temp. Uniformity	+/- 0.5°C at 50°C
Temp. Controller	P.I.D Controller or programmable controller
Temp. Regulator	SSR Type
Capacity (L)	48.5 litres
Circulating bath (optional)	Circulating apparatus
Soxhlet bath (optional)	Ring hole 8ea. Clamp holder
Shaking bath (optional)	Shaking Apparatus
Shaking Range	35-300rpm
Shaking Control	Feed-back control system
Shaking Display and Setting	Digital setting
Safety Device	Electrical Leakage Breaker and Overheating Protection
Material Interior	Stainless steel
Material Exterior	Powder coated steel

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au