

PERFORMER LABORATORY CHEST FREEZERS (-25°C to -15°C)

LABEC supplies 2 models in the **Performer Range of Laboratory Chest Freezers** ranging from 300 litres to 675 litres featuring flat top stainless steel lids and digital display, these freezers provide fast and efficient freezing capacity.

The Labec freezers are designed for general purpose storage of non-critical products where alarms and temperature monitoring are not generally required. Lockable swivel castors fitted to the base make moving the freezer into place easy and general cleaning more accessible.

Controlled temperature conditions within all models using a digital temperature controller with Hi/Lo temperature memory and lift up handle with integrated lock and keys as standard.

Durable PVC coated storage baskets included and high load limit increases the storage capacity.



300 Litre Model



675 Litre Model

Serving the Scientific, Medical and Research industries since 1945

LABEC

LABEC

Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au

SPECIFICATIONS

PERFORMER LABORATORY CHEST FREEZERS (-25°C to -15°C)

Model	PLCF300	PLCF675
Capacity (L)	300	675
Ext. Dim. (mm) W x D x H	1012 x 715 x 903	2055 x 715 x 903
Temp. Controller	Digital Display, Hi/Low Memory	
Temp. Min	-15°C	
Temp. Max	-25°C	
Temp. Setting Accuracy	+/-0.1°C	
No of baskets included	2	3
Doors	Flat top stainless steel	
Refrigerant	R134a (CFC Free)	
Options	Hi/Low audible and visual alarm for connection to Building Maintenance System.	
Power (W)	175	425
AMPS	1.2A	2.9A
Weight (kg)	65	107

Serving the Scientific, Medical and Research industries since 1945

LABEC

LABEC

Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au