



Labec supplies Laboratory Freezers ranging from 300 litres to 1400 litres in one and two door configurations, with solid or double glazed glass doors for easy visibility and convenience. The Labec laboratory freezers are designed for general purpose storage of non-critical products where alarms and temperature monitoring are not generally required. Lockable swivel castors fitted to the base make moving the freezer into place easy and general cleaning more accessible. All models feature self closing doors (with internal lighting). Adjustable open wire white plastic coated full width shelves. Controlled temperature conditions within all models using a digital temperature controller with Hi/Lo temperature memory and BMS alarm contacts as standard.

Specifications

Model	PLFGD300	PLFSD300	PLFSD600	PLFSD900	PLFSD1400
Capacity (litres)	300		600	900	1400
Int. Dim. (cm) HxWxD	150x45x40		122x68x66	122x100x66	155x124x65
Ext. Dim. (cm) HxWxD	184x60x64		200x76x76	200x110x66	215x143x86
Temperature Controller	Digital Temperature Controller				
Temp. Min	-10°C				
Temp. Max	-24°C				
Refrigerant	134a CFC Free				
Number of Doors	1			2	3
Door Material	Glass	Solid	Solid	Solid	Solid
Standard Shelving	6 x fixed shelves		5 x adjustable shelves		15 adjustable shelves
Castors	Standard				
Options	Additional shelves, BMS monitoring, temperature monitoring and logging, glass doors on some models.				
Power Supply	10AMP/240V/50Hz				
Weight (kg)	85		160	210	195