



**Laboratory Equipment Pty Ltd**

## PERFORMER STAINLESS REFRIGERTOR



Manufactured using high grade 304 stainless steel. Interior, exterior and shelving, providing a professional and hygienic environment. Adjustable shelves, split door system, suitable for an assortment of situations. With a top mounted refrigeration system to provide maximum internal storage. Fitted with an electronic thermostat with easy read digital display and powered from single phase 230-240V power sources allowing for hassle free installation.

### Specifications

Model	PLRSS-560	PLRSS-1200
Int. Dim. (cm) HxWxD	195x70x48	195x122x60
Ext. Dim. (cm) HxWxD	205x68x80	205x132x80
Temp. Min	+2°C	
Temp. Max	+10°C	
Capacity (L)	560	1217
Number of Shelves	4	8
Power (W)	2400	
Material Interior and Exterior	Stainless Steel	
Temperature Controller	Digital	
Weight (kg)	120	160

For Further Information contact Laboratory Equipment on:

Email: [sales@labec.com.au](mailto:sales@labec.com.au)

Phone: 02 9560 2811

Fax: 02 9560 6131

[www.labec.com.au](http://www.labec.com.au)