



### Features

- Compliant with Quality Care Pharmacy program (QCCP) Cold Chain Standards
- External Digital temperature controller
- Audible/Visual/BMS signals
- Hi/Lo Temperature memory, Hi/Lo Temperature alarm
- External alarm facility for connection to building maintenance system
- Adjustable/removable shelves, Internal light , Key lock door

### Specifications

Model	PVR-130	PVR-374	PVR440*
Capacity	130	374	440
Int. Dim. (cm) WxDxH	44x44x70	53x45x150	53x50x170
Ext. Dim. (cm) WxDxH	54x54x88	60x65x193	68x65x210
Temp. Controller	Digital Controller		
Temp. Min	+2°C		
Temp. Max	+8°C		
Temp. Setting Accuracy	+/- 0.1°C		
Temp Fluctuation	+/- 2°C		
Temp. Sensor	LAE FM4		
Refrigerant	CFC Free		
Interior Material	Easy Clean Plastic Formed Lining		
Exterior Material	Powder Coated Steel		
Door Material	Solid	Solid or Glass	Glass
Shelves (adjustable)	3	4	5
Options	cable entry port, extra shelves, low temperature safety thermostat		
Power (W)	200	220	300
Power Supply	240V / 50Hz		

\* PVR440 pending approval with QCCP



Approved by the Quality Care Pharmacy Program and The Pharmacy Guild of Australia.