



Laboratory Equipment Pty Ltd

REFRIGERATED WATER JACKETED INCUBATORS



Features

- Stainless Steel Interior
- Stainless Steel Exterior
- Temperature range - 0°C to +60°C
- Temperature uniformity +/-0.2°C
- cooled water in the jacket is pumped through the jacket to maintain uniformity
- Digital microprocessor based PID temperature controller
- Options include cycling temperature and lighting, timer, safety thermostat, extra shelves, fan forced for increased temperature uniformity and rapid temperature recovery

Specifications

Model	IWJC18	IWJC24	IWJC36
Int. Dim. (cm) HxWxD	46x38x38	60x50x50	90x60x50
Ext. Dim. (cm) HxWxD	80x60x60	130x60x68	160x88x88
Temp. Min	0°C	0°C	0°C
Temp. Max	60°C	60°C	60°C
Capacity (L)	66	160	300
Power (W)	400	400	400
Weight (kg)	70	90	110

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au