



**Laboratory Equipment Pty Ltd**

## SAND BATH



### Features

- Temperature to +300°C
- All stainless steel construction
- Simmerstat controller
- Removable stainless steel sand tray

### Specifications

Model	SB1	J-HPS1	J-HPS2
Tray Dim. (cm) HxWxD	45x30x4	45x30x5	60x30x5
Temp. Min	Ambient+5°C		
Temp. Max	300°C		
Power (W)	2400	3000	

For Further Information contact Laboratory Equipment on:

Email: [sales@labec.com.au](mailto:sales@labec.com.au)

Phone: 02 9560 2811

Fax: 02 9560 6131

[www.labec.com.au](http://www.labec.com.au)