



Laboratory Equipment Pty Ltd

Shaking Incubator – Small & Medium Size



A range of small and medium sized bench top shaking incubators. Featuring transparent viewing window for easy visibility which minimises the need to open the chamber.

Specifications

Model	M-TU-400	M-TU-400S	M-TU-453	M-TU-454
System	Forced Air Circulation			
Shaking System	Orbital			
Temp. Min	Ambient+5°C			
Temp. Max	70°C			
Temp Accuracy	+/- 0.1°C			
Temp Uniformity	+/- 1°C (at 37°C)			
Temp Controller	PID			
Rate of Shaking	20-250rpm	30-400rpm	30-300rpm	30-400rpm
Shaking Width	25mm			
Internal Dimensions (WxDxH)(mm)	450x380x260	450x380x260	510x40x430	510x540x430
External Dimensions (WxDxH)(mm)	460x490x430	460x490x430	520x650x500	520x650x500

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au