



Laboratory Equipment Pty Ltd

SHAKING INCUBATORS –Economy Benchtop



Economy range of bench top shaking incubators, heat only and heat and refrigerated models available.

Features

- Powerful triple eccentric drive mechanism
- PID Microprocessor control
- Sealed control panel
- Large easy to read LCD display
- Electronic timer from 1 minute to 500 hours
- Non - volatile memory function
- Stainless steel interior

Specifications

| Model | ZHWY-200B | ZHWY-100B | ZHWY-100C | ZHWY-100D |
|-------------------------|---|---|------------------------|---|
| Performance | Heated & Refrigerated | | Heated Only | |
| Platform Size (mm) | 450x400 | 400x370 | 280x220 | 340x370 |
| Temp. Min | 4°C | | 30°C | |
| Temp. Max | 60°C | | 60°C | |
| Temp Accuracy | +/- 0.1°C | | | |
| Temp Uniformity | +/-1°C | | | |
| Shaking Speed (+/-1rpm) | 30-400rpm | | 30-300rpm | 30-600rpm |
| Timer | 1 minute to 500 hours | | | |
| Capacity | 50ml x 5 100ml x 5 250ml x 4 500ml x 3 | 50ml x 4 100ml x 4 250ml x 3 500ml x 3 | 50ml x 12 100ml x 9 | 50ml x 4 100ml x 4 250ml x 3 500ml x 3 |
| Power (W) | 581 | 566 | 307 | 566 |
| Weight (kg) | 132 | 94 | 54 | 94 |

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au