



Laboratory Equipment Pty Ltd

SHAKING INCUBATORS –Heat Only



Specifications

Model	J-SISN-REC	J-SISN- ROT
Shaking Type	reciprocating action	rotation action
Chamber Size (mm)WxDxH	500x500x290	
Ext. Dim. (mm) WxDxH	540x615x460	
Platform size(mm)	400x400	
Temp Range	Ambient+5°C - 60°C	
Temperature Controller	PID Controller	
Capacity	250-300 ml erlenymer flasks x 23 ea /500ml erlenymer flasks x 16 each	
Temperature Accuracy	+/- 0.1°C at 37°C	
Temperature Uniformity	+/- 0.5°C at 37°C	
Temperature Regulator	SSR Type	
Thermocouple	Pt 100 ohms	
Shaking speed	35-300 rpms	
Material interior	stainless steel plate	
Material exterior	steel plate with powder coated finish	
Weight (kg)	35.2kgs	

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au