

SHAKING INCUBATORS – Large Capacity



Model	J-SCR	J-SCI
Shaking Type	Rotation or reciprocating type	
Environmental Temp range	+5°C to 25°C	
Ext. Dim. (cm) WxDxH	1,130x810/60x950	1,060 x 850 x 800
Temp. range	-10°C to C	Ambient +5°C to 60°C
Temp. accuracy	+/- 0.5°C at 37°C	
Temp. uniformity	+/- 1.0°C at 37°C	
Temp controller	P.I.D Controller	
Temp regulator	SSR type	
Thermocouple	Pt 100 ohms	
Heater (W)	1000	800
Shaking Controller	Digital setting & display	
Shaking speed (rpm)	35-300	
Shaking capacity (load)	250ml Erlenmeyer flask x 42 each or 500ml Erlenmeyer flask x 30each (choice)	
Motor Hp	1/16	
Refrigeration Compressor Hp	3/4	NO
Circulation Fan and Motor	diameter x 100mm x 2ea. 55W	
Insulation	Non flammable polystyrene foam	
Observation Window	Reinforced glass 5T	
Lid	Steel plate with powder coated finish	
Safety Device	Electric Leaker Breakage (E.L.B) overheating protection system	
Shelves	NO	NO

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au