



Laboratory Equipment Pty Ltd

SHAKING WATER BATH



Specifications

Model	J-SWB35	J-SWB60
Shaking Type	Reciprocating	
Overall size (mm) (WxDxH)	560x410x380	700x 480 x 380
Internal Dimensions of bath (WxDxH)	400x350x250	600x400x250
Operating temperature	+5°C to 35°C	
Temp. Range	Ambient+5°C to 80°C	
Temp. Accuracy	+/- 0.1°C at 50°C	
Temp Uniformity	+/- 0.5°C at 50°C	
Temp Controller	P.I.D setting/digital display	
Temp Regulator	SSR Type	
Thermocouple	Pt 100ohms	
Shaking Speed range	30-200rpm	
Shaking stroke	30mm	
Shaking capacity	250ml x 9each	500ml x 12 each
Shaking plate (mm)	320 x 310	520 x 360
Heating time	Reaches 50°C in 30 mins	
Display	Digital	
Capacity (L)	35	60
Heater (sheath)	3000W	
Motor	1/8 HP	
Interior material	Stainless Steel	
Exterior material	Powder Coated Steel	
Safety Device	Electric Leaker Breakage (ELB)	
Power (W)	600	500

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au