



**Laboratory Equipment Pty Ltd**

## SIEVE SHAKER



### Specifications

<b>Model</b>	<b>J-893</b>
Shaking Type	Ro-tap type
External Dimension (mm) WxDxH	680x530x670
Ro-tap speed	240 rpm
Hammer speed (up/down)	130 rpm
Motor Hp	1/3 Hp
Sieve Quantity	8 inch diameter/2 inches high x 6 each
Sieve size	8 inch diameter/2 inches high
Timer	60 minutes
Weight (kg)	75.2

For Further Information contact Laboratory Equipment on:

Email: [sales@labec.com.au](mailto:sales@labec.com.au)

Phone: 02 9560 2811

Fax: 02 9560 6131

[www.labec.com.au](http://www.labec.com.au)