



Laboratory Equipment Pty Ltd

SOXHLET BATHS



Features

- All stainless steel construction
- 3 or 6 hole configuration of 100mm concentric rings - 5 rings per hole.
- Removable concentric rings
- Aluminium supports arms at rear
- Simmerstat control for temperature
- Automatic overheat/water low cut off protection requiring user to restart after fault.

Specifications

Model	SX3	SX6
Int. Dim. (cm) HxWxD	40x17x10	80x17x10
Ext. Dim. (cm) HxWxD	43x26x23	86x26x23
Temp. Min	Ambient+5°C	Ambient+5°C
Temp. Max	100°C	100°C
Capacity (L)	7	15
Power (W)	1200	2000
Weight (kg)	8	11

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au