



**Laboratory Equipment Pty Ltd**

## **SPLIT TUBE FURNACES**



### **Features:**

- Stainless steel construction
- Designed for to suit applications such as mechanical testing apparatus
- Up to 1300°C
- Light weight ceramic fibre insulation for rapid thermal cycling
- Remote control units
- Digital microprocessor based PID temperature controller
- All units custom made to suit your requirements
- Optional multi stage fully programmable controllers available

1300°C models available on request

\*\* Vertical tube furnaces made to order

\*\*\* Split tube furnaces made to order

Specifications are subject to change without notice

Custom built sizes available on request

For furnaces with operating temperatures above 1300°C please contact our office

For Further Information contact Laboratory Equipment on:

Email: [sales@labec.com.au](mailto:sales@labec.com.au)

Phone: 02 9560 2811

Fax: 02 9560 6131

[www.labec.com.au](http://www.labec.com.au)