



Laboratory Equipment Pty Ltd

STEAM STERILISER



Features

- Stainless steel vessel
- Simmerstat control
- Gabled lid with vent
- Constant water level from mains supply
- Safety cut-out device
- Rapid steam production

Specifications

Model	SS1	SS2
Int. Dim. (cm) WxLxD	40x40x50	40x20x20
Ext. Dim. (cm) HxWxD	90x55x60	25x50x25
Temp. Min	Ambient+5°C	
Temp. Max	100°C	
Capacity (L)	80	16
Power (W)	3600	2400
Weight (kg)	30	22

For Further Information contact Laboratory Equipment on: Email: sales@labec.com.au

Phone: 02 9560 2811 Fax: 02 9560 6131 www.labec.com.au