



Laboratory Equipment Pty Ltd

TISSUE FLOAT BATHS



Features

Manufactured with aluminium casting, the Tissue float bath is used for removing substances such as paraffin's wrinkle and foam (bubble) in pathology and hospital laboratories.

Specifications

Model	J-NBT
Size Dia x H(mm)	220x60
Temp. Min	Ambient+5°C
Temp. Max	60°C
Heater(sheath) W	300
Weight (kg)	2.75

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au