



Laboratory Equipment Pty Ltd

ULTRA LOW TEMPERATURE DEEP FREEZER CHEST MODEL -40°C to - 86°C



Labec now supply a range of low temperature freezers. Our Chest Model Ultra Low Temperature Deep Freezer is practical and convenient and designed to make good use of the laboratory area.

Features

- Saving Power Consumption
- Optimum interior design allows stable temperature.
- Urethane injection allows longer durability.
- RS-232: You can check/store operating condition of deep freezer and temperature can be controlled on the remote site
- Is equipped with Automatic Voltage Boosting System so Deep Freezer is protected at all times.

Options

- Chart Recorder

Specifications

Model	I-DF9007	I-DF9010	I-DF9014	I-DF9017	I-DF9020
Int. Dim. (mm) HxWxD	670x570x500	670x870x500	670x1170x500	670x1470x500	670x1770x500
Ext. Dim. (mm) HxWxD	1030x1280x790	1030x1590x790	1030x1880x790	1030x2180x790	1030x2480x790
Refrigeration	2 x 1/2 HP				2 x 1 HP
Temp. Min	-86°C				
Temp. Max	-40°C				
Capacity (L)	191	295	392	492	593
Weight (kg)	211	236	261	286	311

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au