



Laboratory Equipment Pty Ltd

PROGRAMMABLE UNIVERSAL CENTRIFUGE



Features

- LED display speed and time
- 10 memory programs for presetting parameters
- Microprocessor controller. Auto accelerating and decelerating system
- Auto cover lid lock. Auto brake system for presetting
- Safety device for rotor imbalance
- Auto buzzer system when operation is finished
- Max speed 6000 rpms. Max RCF 3980x g
- Power cut off when the cover is open. Auto balance device
- Air auto cooling circulation.ABS outer covering; see through window
- Other rotors and accessories available
- Dimensions WxDxH (mm) 290x320x245

Specifications

Rotor	Time & Speed	Max RCF (x g)	Capacity (ml x pcs)	Tube Size Max (dia x H)	Tube Size Min (dia x H)
RA-1508	0-99 min.	3,940	15 x 8	16.7 x 122	15 x 90
RA-1521S		3,940	15 x 12	16.7 x 122	15 x 90
RA-5004	500-6000 rpm	3,980	50 x 4	29.5 x 117	27.5 x 90
RA-5006		3,980	50 x 6	29.5 x 117	27.5 x 90
RA-1520		3,140	1.5 x 20	11 x 48	10 x 40

For Further Information contact Laboratory Equipment on:

Email: sales@labec.com.au

Phone: 02 9560 2811

Fax: 02 9560 6131

www.labec.com.au