

DIGITAL UNIVERSAL CENTRIFUGE

LABEC offer a compact bench top universal centrifuge with the capacity for rotors ranging from 1.5ml to 50ml and speeds of 500-6000rpm's.

Features:

- *LED display speed and time
- *10 memory programs for presetting parameters
- *Microprocessor controller
- *Auto accelerating and decelerating system
- *Auto cover lid lock
- *Auto brake system for presetting
- Safety device for rotor imbalance
- *Auto buzzer system when operation is finished
- *Max speed 6000 rpms
- *Max RCF 3980x g
- *Power cut off when the cover is open



Rotor	Speed	Max RCF	Capacity	Tube Size	Tube Size
	rpm's	(x g)	(ml x pcs)	Max (dia x H)	Min (dia x H)
RA-1508	50-8000rpm	3,940	15 x 8	16.7 x 122	15 x 90
RA-1521S	50-8000rpm	3,940	15 x 12	16.7 x 122	15 x 90
RA-5004	50-8000rpm	3,980	50 x 4	29.5 x 117	27.5 x 90
RA-5008	50-8000rpm	3,980	50 x 8	29.5 x 117	27.5 x 90
RA-1520	50-8000rpm	3,140	1.5 x 20	11 x 48	10 x 40

Serving the Scientific, Medical and Research industries since 1945

LABEC



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd
email: sales@labec.com.au
Ph: 02 9560 2811 • Fax: 02 9560 6131
www.labec.com.au