

WATER DISTILLER

LABEC Water Distiller uses a cool water supply to convert tap water into high grade distilled water for use in all general applications where distilled water is required.

The cold water can be re used from an optional storage tank pumped back to the distiller. The distiller uses minimal space and can run continuously providing a constant supply of distilled water. Set up and installation is extremely easy and does not require any fixed plumbing or permanent fixture.

The interior and exterior are made from anti corrosive stainless steel, a water leveller automatically cuts off power supply to stop operation when the water supply is low.



SPECIFICATIONS

Model	J-WD	J-WD1
Environmental temp. range	+5°C to +40°C	
Heater (sheath)	3000	6000
Output per hour (litres)	4 litres	8 litres
Power (W)	15AMPS/240V	25AMPS/240V
Material	Stainless Steel	
Weight (kg)	6.5kg	11.7kgs

Serving the Scientific, Medical and Research industries since 1945

LABEC



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au