



**Laboratory Equipment Pty Ltd**

## WATER JACKETED INCUBATORS



### Features

- Designed for applications where non-drying of samples is essential
- Stainless Steel Interior and Exterior
- Temperature range - ambient +5°C to +70°C
- Digital microprocessor based PID temperature controller
- Stainless Steel interior water jacket
- Minimal temperature fluctuation +/-0.1°C
- Options include fan forced air circulation, extra shelves, timer and safety thermostat

### Specifications

Model	IWJ14	IWJ18	IWJ24	IWJ36	IWJ48
Int. Dim. (cm) HxWxD	35x30x30	46x38x38	60x50x50	90x60x50	120x60x50
Ext. Dim. (cm) HxWxD	65x50x40	80x58x48	97x70x60	120x80x60	150x80x60
Temp. Min	Ambient+5°C	Ambient+5°C	Ambient+5°C	Ambient+5°C	Ambient+5°C
Temp. Max	70°C	70°C	70°C	70°C	70°C
Capacity (L)	32	66	160	300	375
Power (W)	150	150	150	400	400
Weight (kg)	65	90	120	120	130

For Further Information contact Laboratory Equipment on:

Email: [sales@labec.com.au](mailto:sales@labec.com.au)

Phone: 02 9560 2811

Fax: 02 9560 6131

[www.labec.com.au](http://www.labec.com.au)